SomaSuite™ v1.0.x

Download and Installation Instructions

Somalogic

Overview and Download Instructions

Overview

SomaSuite was developed by SomaLogic in partnership with NEC to assist in exploratory visualization and analysis of proteomic data produced by the SomaScan Assay. The software has five sequential functions: 1) load data files, 2) define sample groups, 3) perform statistical analysis, 4) create plots based on sample group and statistical analysis filters, and 5) format results. SomaSuite is a standalone software package that is intended for installation under Windows XP or Windows 7 and is an English language only format.

Download Instructions

These instructions are for SomaLogic clients who may not regularly use FTP to store and retrieve files. If you have a preferred FTP client, the necessary information is listed below and you may wish to use your client to download the software and proceed to the **Installation Instructions**.

FTP: <u>ftp.somalogic.com</u>

User: soma1234 PW: Somalogic1!

Secure FTP is a safe and effective method for transferring large files from host to host using an established File Transfer Protocol over the Internet. Directly opening the link above will typically initiate the FTP using your default Internet browser which may result in an unstable connection.

Step Download Instructions

Preview:

1 Open Windows Explorer

You may use Start, Documents or Start, My Documents



2 Enter the ftp address in the address bar:

ftp.somalogic.com will resolve to: ftp://ftp.somalogic.com



Overview and Download Instructions

Download Instructions Step


Preview:

3 Enter the User name and Password: User name: soma1234

Password: Somalogic1!




4 Copy the SOMASUITE Software folder:



Paste the folder locally: 5 (Desktop, My Documents, or alternative download folder)

> This process typically takes 4-6 minutes depending on download speed.



Overview and Download Instructions


Step Download Instructions

downloaded files

6

Open the local folder to view the

Preview:




Step Installation Instructions

1 Double-click on the install file: SOMASUITE v1.0.x


Select Next

Preview:



2 The Welcome screen provides an early opportunity to cancel.

Select Next.



Step Installation Instructions

The SOMASUITE End User License contains important information regarding the terms of use of the Software. Please read the agreement.

Select I Agree

Select Next


Preview:



4 Optional: Select Browse to change the default installation path

Optional: Select Disk Cost to calculate the available space and space used by the software installation.

Select Next




Step Installation Instructions

5 The confirm screen provides a final opportunity to cancel the installation.


Select Next to install

Preview:



6 The installation typically takes 2-3 minutes. When installation is complete:

Select Next




Step Installation Instructions

7 Please contact SomaLogic if the installation fails to Complete.
Otherwise:

Select Close

Preview:




8 A Complete installation will result in a SOMASUITE Icon on the Desktop.

Double-click on the SOMASUITE Icon to start the application.



9 Select Getting Started

Review the brief SOMASUITE tutorial to familiarize yourself with the basics of the software. Additional details and instructions can be found in the User's Manual.



Step Installation Instructions

To load data and begin analysis:Select Start New Analysis

Preview:

